



The aLIGO Particulate Evaluation Tool (PET)

AUTHOR(S)	DATE	Document Change Notice or Approval
Kate Gushwa and Calum Torrie	January 10, 2014	See LIGO DCC Record

Table of Contents

- 1 INTRODUCTION 3**
- 2 SCOPE 3**
- 3 IMPLEMENTATION..... 3**
 - 3.1 Introduction Period 3
 - 3.2 Pilot Trial..... 4
 - 3.3 Team Trial..... 4
 - 3.4 Full Adoption 4
- 4 SURFACE PARTICULATE COLLECTION..... 5**
 - 4.1 Required Parts and Tools 5
 - 4.2 Cleaning 6
 - 4.3 Sample Collection Method 7
 - 4.4 Control Samples..... 9
- 5 STEREO MICROSCOPY & IMAGE ANALYSIS SYSTEM SETUP 10**
 - 5.1 Required Parts and Tools 10
 - 5.2 Cleaning 11
 - 5.3 Microscope Setup..... 11
 - 5.3.1 LED Stand..... 11
 - 5.3.2 Microscope Body 12
 - 5.3.3 X-Y Translation Stage 13
 - 5.3.4 Rough Alignment 14
 - 5.3.5 Acetate Template 15
 - 5.3.6 Total Magnification 16
 - 5.3.7 Microscope Settings Summary 16
 - 5.4 Software Setup 17
 - 5.4.1 ToupView Preferences 17
 - 5.4.2 Camera Settings 17
 - 5.4.3 AmScope Fine Alignment..... 19
 - 5.4.4 Calibration..... 20
- 6 SAMPLE ANALYSIS PROCEDURE..... 22**
 - 6.1 Required Parts and Tools 22
 - 6.2 Sample Preparation and Setup..... 22
 - 6.3 Stitching..... 24
 - 6.3.1 Basic Method 24
 - 6.3.2 Advanced Method 26
 - 6.4 Image Processing 27
 - 6.5 Image Size..... 30



The aLIGO Particulate Evaluation Tool (PET)

6.6	Particle Count	31
6.7	Summary of Steps and Settings	33
7	PARTICLE CLEANLINESS LEVEL CALCULATOR	35
7.1	Raw Data Entry & Results	35
7.2	Definition of Cleanliness Levels	37
7.3	Calculations	38
7.3.1	Areas	38
7.3.2	Binning	39
7.3.3	PCL	39
8	SUMMARY OF SETTINGS	40
8.1	Microscope	40
8.2	ToupView	40
9	TROUBLESHOOTING	42
9.1	Black Camera View	42
9.2	Focusing	42
9.3	Incomplete Stitching	43
9.4	Failed Stitch	44
9.5	Count Plugin Will Not Open	44
9.6	Particles Lost/Gained	44
9.7	Toolbars Disappeared	44

Special Note: The aLIGO Particulate Evaluation Tool is based on a similar technique, called the Particle Cleanliness Validation System (PCVS), developed by staff at the National Ignition Facility (NIF). A special thanks to NIF staff for sharing reports and key information on this system during the visits and discussions with the authors of this document.



The aLIGO Particulate Evaluation Tool (PET)

1 INTRODUCTION

The Advanced LIGO interferometers are extremely sensitive to optical scattering and absorption losses induced by both particulate and hydrocarbon contamination. Because optical performance is key to the success of the interferometers, reducing contamination in the vacuum system is critical.

Cleanrooms are defined by a cleanliness class, which refers to the permitted quantity and size of particles *per volume of air*. Similarly, surfaces in cleanrooms are defined by a cleanliness level, a standard set of particle size distributions criteria *per area*. It is already standard practice to check air quality in LIGO's ISO 5/Class 100 cleanrooms with a handheld particle counter. Although IEST-STD-CC1246 defines surface cleanliness, called Particulate Cleanliness Levels (PCLs), it does not specify exactly how to measure them. Following a consultation with the National Ignition Facility (NIF) on cleanliness practices, a new Particulate Evaluation Tool (PET) has been developed to quantify contamination on mechanical surfaces. The PET has three main components:

- Dry swipe technique and particulate collection tool to take samples from surfaces,
- Sample analysis with a stereomicroscope and ToupView software to count particle areas,
- Excel spreadsheet to compute the particle size distribution and PCL.

This document provides a detailed explanation of how to set up and use the PET, which can be thought of as a white glove test or the surface equivalent of a particle counter for cleanroom air. It should be noted that this technique is NOT meant to be used for optical surfaces.

2 SCOPE

The current method for measuring PCLs is to place witness wafers and optics in chamber, and then perform a particle size distribution analysis on the contamination that is collected. An automated dark field microscope provides precise particle counts for major installation work and pump/down vent cycles. The major disadvantage of this approach is the incredibly long lead time. Wafers and optics spend weeks to months in chamber, days being shipped from the observatories, and days to weeks at Caltech before analysis.

The Particulate Evaluation Tool is the perfect complement to the witness wafers and optics. The entire process can be quickly performed on site so install teams can immediately assess the condition of the area they are working in. Long-term, the PET will be used to track changes in contamination during different installation activities, determine the effectiveness of cleaning practices, verify chambers have been properly cleaned, and measure overall progress towards reaching LIGO's cleanliness goals – PCL 65 for vertical optical surfaces and PCL 100 for horizontal mechanical surfaces. Refer to [LIGO-T1300511](#) for additional information on aLIGO particulate contamination requirements.

3 IMPLEMENTATION

A phased implementation strategy will be used to effectively incorporate the new PET system into current aLIGO activities.

3.1 Introduction Period

The first step was to set up the microscope in a dedicated area, and calibrate its software for sample analysis. Members of the Contamination Control Working Group (CCWG) from Caltech were on site in November 2013 to assist with the initial set up and to train at least one person to use the system. Travis Sadecki became the Hanford PET analyst for the



The aLIGO Particulate Evaluation Tool (PET)

CCWG. Bryan Smith and Stuart Aston will share the responsibility at Livingston. Future users will be introduced to the PET during presentations on contamination control and morning meetings.

3.2 Pilot Trial

From November 2013 to early January 2014, the newly trained analysts will swipe a variety of clean and dirty surfaces from multiple locations all over the sites. A few of these samples should be taken of the same surface before and after specific activities, such as cleaning or assembly work. The pilot trial will be used to assess the impact (if any) of sample collection on other work activities, and the actual time required to both take and analyze samples. In theory, samples can be taken in parallel to ongoing work with minimal interference, even in chamber. The pilot trial will highlight any errors or shortcomings in the procedure, allowing ample time for improvement before widespread implementation of the PET at both observatories.

During this period, the analysts will become familiarized with the PET procedure and begin to recognize patterns in the initial data, such as the effects of specific activities or areas in need of additional cleaning. The analysts will present their findings to the CCWG.

3.3 Team Trial

In January 2014, the CCWG and PET analysts will teach select members of the installation teams how to sample surfaces with the swipe tool. The trial teams will regularly collect samples before AND after events requiring cleaning per the “clean-as-you-go” procedure, [LIGO-G1300777](#):

- Assembly (at key steps in the assembly process)
- End assembly (before storage/transportation)
- Pre-Installation (after removing from storage/transportation)
- Post-Installation
- Pre-Chamber close

Though sample collection will become a shared task, sample analysis and PCL tracking will continue to be the responsibility of the site’s PET analyst. The analysts will report on surface cleanliness during morning meetings and in aLOGs.

3.4 Full Adoption

If the team trials are successful, the PET will be fully adopted at both observatories. The PET analysts will lead the effort, and advise when and where samples should be collected. Anyone who works in cleanrooms, including cleaners and visitors, could be expected to collect swipe samples.

The aLIGO Particulate Evaluation Tool (PET)

4 SURFACE PARTICULATE COLLECTION

This section describes the tools and method used to collect contamination samples from Class A and Class B surfaces. *The PET should NEVER be used on wet surfaces or optical surfaces.*

4.1 Required Parts and Tools

The Particulate Collection Tool, or swipe tool, was developed specifically for taking dust samples from tables and instruments inside cleanrooms and chambers. In addition to a swipe tool and foam cutting guide plate, Systems provided both observatories with an initial set of white membrane filters, filter forceps, PetriSlides, and sponges. A few general lab supplies, such as a ruler, are also needed for sample collection. The full list of required parts and information for ordering additional supplies is given in Table 1.

Keep the surface particulate collection tools in a dedicated box or area so that specialty items, like the filter forceps, do not get mixed in with general lab supplies.




Figure 1. Swipe tool and accessories

Table 1. Parts required for surface particulate collection

Part	P/N	Vendor	Purchasing Responsibility
Swipe tool	LIGO-D1300697	N/A	Contact authors for additional tools
Foam cutting guide plate	LIGO-D1301010	N/A	Contact authors for additional plates
Membrane filters	ATTP04700 (100/pack)	Millipore , Amazon	Sites to re-order as required
PetriSlides	PD1504700 (100/pack)	Millipore , Amazon	Sites to re-order as required
Filter forceps	XX6200006P (3/pack)	Millipore , Amazon	Sites to re-order as required
UltraSOLV sponge	HT4754	VWR	Sites to re-order as required
X-acto knife	XZ3601	Amazon	Sites to re-order as required
Replacement blades	XZ211	Amazon	Sites to re-order as required
Cutting mat	GBM0812	Amazon	Sites to re-order as required
Flashlight array	LIGO-D1300223	N/A	Contact authors for additional arrays
Vectra Alpha 10 wipes	TX1010	VWR	Sites (general lab supplies)
IPA	N/A	N/A	Sites (general lab supplies)
Metal ruler	N/A	N/A	Sites (general lab supplies)

The aLIGO Particulate Evaluation Tool (PET)

4.2 Cleaning

The swipe tool was sent to the observatories already assembled, and should not be completely disassembled for cleaning. Remove the sponge piece from the roller section, and thoroughly scrub the entire swipe tool with an IPA soaked Vectra Alpha 10 wipe. When the swipe tool is dry, insert the sponge again.

An IPA soaked wipe is also used to clean the PetriSlides, cutting mat, and the outside of the membrane filter container. If particulates are found inside the filter container, gently wipe the inside of the plastic lid and around the base, taking care not to wet the filters. The filters should not be removed to clean the container. For cleanliness, each membrane filter is separated by a piece of blue paper.

The filter forceps, X-acto knife, and foam cutting guide plate should be Class B cleaned.

The sponge will become dirty or damaged over time, and will need to be replaced. To cut a new piece, lay a clean cutting mat on a work table. Place the foam block on the mat, and the guide plate on top of the foam. Align the edges of the guide plate to the corners of the foam. Cut the foam with the X-acto knife in the order indicated in Figure 3. If the cut sponge piece is not completely separated from the foam block, then remove the guide plate and sever the connection. It is a good practice to cut multiple sponges at a time for later use.




Figure 2. Swipe tool parts (ISO exploded view)




Figure 3. Foam cutting guide plate (top view)




Figure 4. Cutting sponge with X-acto knife and guide plate

The aLIGO Particulate Evaluation Tool (PET)

4.3 Sample Collection Method

The following is a step-by-step explanation of the surface particulate collection method. The technique should be practiced at least once before taking samples for particle counting. New users are also encouraged to tear a filter membrane to get a better feel for how fragile the filters are.

- 1) Using clean forceps, remove and discard the top blue separation paper from the stack of membrane filters.

Note: If a white filter is at the top of the stack, it may be contaminated and should be discarded.

- 2) Remove a clean membrane filter by gently picking up a small portion near the edge with the forceps, and visually inspect it for contamination. If any particles are visible on the filter, it should be discarded.

Note: Contaminated filters will yield inaccurate PCLs. To avoid contaminating the filter, minimize the time it spends outside of its original container and PetriSlide.

- 3) Press down on the swipe tool's clamp, insert the filter between the tool's body and clamp, then release the clamp to secure the filter in place. The widest part of the filter should be aligned with the sponge, as denoted by the red line in Figure 6.

Particles are collected on the portion of the filter directly over the sponge, which is referred to as the sample collection area.

Note: The forceps can leave a mark or tear in the filter, which may be miscounted as a large particle. To avoid inaccurate counts, do not touch the sample collection area.

- 4) If possible, use a flashlight array to illuminate the surface that will be sampled.
- 5) Place a ruler on or near the site that will be sampled to measure the length of the swipe. The swipe distance must be at least 12 inches. The longer the swipe, the more accurate the PCL calculation. The same filter can be used to swipe a surface more than once if there is insufficient space to take a single 12 inch long swipe. For




Figure 5. Removing blue separation paper




Figure 6. Inserting membrane filter into swipe tool



Figure 7. Side view of swipe tool and filter

The aLIGO Particulate Evaluation Tool (PET)

example, four shorter swipes can be taken along the base of a suspension. For multiple swipes, add the length of each swipe.

Note: It is crucial for wipe distance to be measured as precisely as possible, not just guesstimated. PCLs will be tracked over time to improve our understanding of what events affect cleanliness, measure the effectiveness of contamination control practices, and evaluate how well teams have cleaned chambers. Accurate measurements are needed to avoid apples to oranges comparisons.

- 6) Use a finger to hold the trailing edge of the filter to the swipe tool body without touching the sample collection area. The correct technique is demonstrated in Figure 9.

Note: The procedure did not originally include this step, as depicted in Figure 10. The result was particles spread over larger areas, which required either significantly more time for image processing or the cropping of large portions of the filter to save time. Constraining the sample collection area solves this issue.

- 7) While holding the swipe tool at a roughly 45 degree angle, apply a light vertical force and slowly drag the filter flat across the surface. Less than 1 pound of force is needed. Applying too much force can scratch, pockmark or tear the filter.

Note: Indentations or tears in the sample collection area will make particle counting difficult or impossible. A filter can still be used for counting if there are fewer than TBD small scratches (less than TBD mm) present anywhere within the sample, or if small rips (less than 1 mm) are found around the edges of the sample collection area. Refer to Section 6 for additional information. If unsure whether or not a sample should be discarded, contact the site PET analyst or one of the authors.




Figure 8. Ruler to measure swipe length




Figure 9. Correct swiping technique




Figure 10. Incorrect swiping technique

The aLIGO Particulate Evaluation Tool (PET)

- 8) Check the particle removal efficiency of the swipe. A flashlight array is useful for this task. The filter should collect about 85% of particulates (though it can become saturated in very dirty areas). Increase the vertical force applied to the tool or swipe over the same area a second time if needed. If large particles or fibers are left on the surface after swiping, add them to the sample with tweezers after Step 10.
- 9) Release the swipe tool clamp, and gently remove the filter from the swipe tool with the forceps. Do not touch the sample collection area.
- 10) Place the filter into a clean PetriSlide, with the sample particulates facing up.
- 11) Label the PetriSlide, including as much information as possible. At the very least, the label should have the date, location, and swipe length.
- 12) Thoroughly wipe down the area that was swiped with an IPA wetted Vectra Alpha 10 wipe.




Figure 11. Sample particles on membrane filter




Figure 12. Filter in PetriSlide

4.4 Control Samples

The following control samples should be prepared at each site:

- Blank (un-swiped membrane filter)
- Class A surface fresh out of the oven or storage
- Class A surface just after it has been cleaned with an IPA wetted Vectra Alpha 10 wipe
- “Dirty” surface, such as an optic table in a cleanroom

The control samples are used in Section 6 to help the PET analysts learn the steps involved in the imaging and counting procedure, and to provide a baseline calibration between users at the same site. The control samples are also an example of the PCLs that can be achieved by cleaning.



The aLIGO Particulate Evaluation Tool (PET)

5 STEREO MICROSCOPY & IMAGE ANALYSIS SYSTEM SETUP

The Livingston and Hanford observatories and Caltech LIGO Laboratory have been provided with an identical AmScope stereo microscope, digital camera, LED stands, Lenovo laptops, and image analysis software. This section summarizes the main parts of the microscope and their purposes, and also describes how to calibrate the hardware and software for the analysis of filter samples. ***This setup should only need to be performed once (if the microscope is not moved).***

5.1 Required Parts and Tools

Collect the items listed in Table 2, and print out the PET Microscope X-Y Stage Template, [LIGO-D1301007](#), on acetate prior to assembling the microscope. When printing the template, make sure to select “actual size” not “fit to page.” Use a guillotine paper cutter to trim the templates along the dotted lines.

Table 2. Parts required for initial microscope setup

Part	P/N	Vendor	Purchasing Responsibility
Microscope + 5MP camera	ZM-1TW3-FOR-5M	AmScope	Contact authors for replacement parts
LED stand	TS130-LED	AmScope	Contact authors for replacement parts
Laptop	T530	Lenovo	Contact authors for replacement parts
Image analysis software	ToupView (64 bit)	ToupTek	Contact authors for replacement CD
Micrometer calibration slide	MR400	AmScope	Sites to re-order as required
X-Y stage	GT100	AmScope	Contact authors for replacement parts
PET custom microscope templates	LIGO-D1301007		
Acetate	3R3028 (100/pack)	Amazon	Sites to re-order as required
Printer	N/A	N/A	Sites (general lab supplies)
Guillotine paper cutter	N/A	N/A	Sites (general lab supplies)
Machined inserts for stand	N/A	N/A	Contact authors for replacement parts
Anti-rotation locks	N/A	N/A	Contact authors for replacement parts
Microscope slides	206A2 (72/pack)	Amazon	Sites to re-order as required
Vectra Alpha 10 wipes	TX1010	VWR	Sites (general lab supplies)
IPA	N/A	N/A	Sites (general lab supplies)
Screwdrivers Phillips #0 Phillips #2 Slotted 3/16"	N/A	N/A	Sites (general lab supplies)
Metal ruler	N/A	N/A	Sites (general lab supplies)
Kapton tape	N/A	N/A	Sites (general lab supplies)

The aLIGO Particulate Evaluation Tool (PET)

5.2 Cleaning

The microscope must be assembled and maintained in a clean environment to prevent contamination of its optics and the samples being analyzed. The alignment and calibration process will have to be repeated if the microscope is moved; therefore, a dedicated area for the AmScope is preferable, such as a flow bench in the optics lab. Select an isolated location with minimal foot traffic because even seemingly minor disturbances to the setup can negatively impact the counting procedure or require re-alignment. Ideally, the space will have enough room for storing the swipe tools and related accessories, AmScope and its unused parts, and growing library of surface samples.

Excluding the lenses, all of the microscope parts arrived on site dirty. All of the microscope parts EXCEPT the lenses need to be wiped down with an IPA soaked wipe. The same applies for the X-Y stage, acetate template, calibration slide, and microscope slides. Take care not to scratch delicate surfaces while cleaning. If the lenses require cleaning, please contact the authors for further instructions.

The ruler, scissors, and screwdrivers required for setup should already be available on site. If these general lab tools are not already Class B cleaned, then wipe them down with IPA. Always keep the microscope covered with aluminum foil or C3 when it is not in use.

5.3 Microscope Setup

5.3.1 LED Stand

The LED stand is a replacement for the stock microscope stand. The new stand has upper and lower sources of illumination, which are controlled by individual knobs on the base of the stand. Mount the upper LED to the arm, and connect its cable to the metal pole, as depicted in Figure 14. Plug the stand in and turn on the power. Confirm that both




Figure 13. Dedicated PET area in Caltech Contamination Lab




Figure 14. LED stand (left view)

The aLIGO Particulate Evaluation Tool (PET)

illuminators are functioning properly. There should be no flickering or burned out LEDs. When analyzing samples, both LEDs should be on the brightest setting. Adjust the position of the upper LED to light the area over the lower LED, and tighten its screws to hold in place.

The arm is the structural element that connects the microscope body to the base. It can move up and down, as well as rotate around the pole. Tightening the clamping knob fixes the position of the arm. The stop unit is a safety feature that prevents the arm from sliding further down the pole. Move the stop unit and arm several inches up the pole, and then lock.

The stage clips and circular glass stage plate are not needed for the PET setup. Remove these parts, and store with the rest of the unused microscope parts.

5.3.2 Microscope Body

The body is the upper part of the microscope that houses the optical components. Loosen the body lock screws on the arm, and insert the microscope body into the ring in the arm, as in Figure 15. Tightening the body lock screws will prevent the body from rotating inside the arm.

Several objective lenses and eyepieces in labelled boxes were provided with the AmScope. Attach the 1X objective lens and 10X eyepieces to the microscope body.

Mount the digital camera on the fixed camera adapter, which has a 10X magnification. The height of the camera is adjusted by rotating the adapter. Position the camera at its lowest position, i.e. as close to the microscope body as possible. Connect the USB cable to the digital camera. For binocular observation, slide the light path selecting lever to “OUT.” Select “IN” to use the digital camera. The eyepieces must be covered with foil to block ambient light when the digital camera is in use.

Anti-rotation locks were created for the




Figure 15. Microscope body with 1 anti-rotation lock (right view)

The aLIGO Particulate Evaluation Tool (PET)

camera by Systems, and are attached to the body using the adapter's locking screws. Only one anti-rotation lock is installed in Figure 15 to provide an unobstructed view of the digital camera and its adapter. In the final configuration, the microscope should have two anti-rotation locks.

In addition to the fixed magnification achieved by the objective lens, the AmScope has a continuously variable magnification range that can be adjusted by turning a knob. The zoom magnification knob is located on the right side of the microscope. Align the 1X marker with the raised line on the microscope body just left of the zoom control knob.

5.3.3 X-Y Translation Stage

As Figure 16 illustrates, the specimen clips on the X-Y translation stage are not only unnecessary but also irksome because they scratch the paint off the top plate. Detach and permanently store the clips. Remove and temporarily set aside the plate stop, its two screws, and the glass plate.

Fit the X-Y stage onto the stand's base with the stage opening directly over the lower LED, as demonstrated in Figure 17. Two custom machined inserts added to the clip holes in the stand's base will hold the X-Y stage in place. One of the inserts has a flat side which must be flush with the side of the stage, otherwise the stage will not be held securely. Push both inserts as far into the holes as possible to avoid blocking movement of the Y-stage. Turn the translation knobs to confirm there are no obstructions to




Figure 16. X-Y Translation Stage




Figure 17. X-Y stage on LED stand (right view)



The aLIGO Particulate Evaluation Tool (PET)

prevent the stage from sliding smoothly in both directions.

The working distance is the distance between the objective lens and the sample, or in this case, the top of the X-Y stage. With the stop unit locked, set the working distance to 4 inches by loosening the clamping knob and moving the arm up or down the main pole.

5.3.4 Rough Alignment

The machined inserts skew the position of the X-Y stage with respect to the stand's base. The digital camera and stage must be aligned to successfully complete the image analysis procedure outlined in Section 6. Since the orientation of the camera is dependent on the configuration of the microscope body and arm, those three components must be properly positioned before rough alignment with the stage.

The centerline of the microscope body should be parallel to that of the arm and the digital camera, as highlighted by the blue dotted lines in Figure 18. If the camera adapter lock screws are tight, the anti-rotation locks will keep the camera parallel to the arm. To adjust the position of the microscope body relative to the arm, loosen the body lock screws and then rotate the body inside the arm ring until the arm and body centerlines coincide.

The microscope body, arm, and camera should now be roughly aligned and considered a single unit. The next step is to align this assembly with the offset X-Y stage. With the stop lock securely in place, the clamping knob can be loosened to allow the arm to swing around the main pole in an arc until the objective lens is positioned directly above the circular stage opening.

Fit the metal ruler into the grooves on the X-Y stage that normally hold the glass stage plate, as shown in Figure 18 and Figure 19. Examine the position of the ruler using the eyepieces. If needed, use the focusing knobs located on the left and right side of the microscope. To change the tension for finer adjustment, hold one of the focusing knobs and turn the other clockwise. Turning counterclockwise will loosen the focusing knob and is not recommended. Moving the stage only along the x-axis, check the alignment of the setup. Make adjustments to the position of the microscope body and arm alignment if necessary. Then, do the same for the y-axis. Fine alignment will be performed in Section 5.4.3.




Figure 18. AmScope alignment (right view)



Figure 19. AmScope alignment (front view)

The aLIGO Particulate Evaluation Tool (PET)

5.3.5 Acetate Template

Store the metal ruler and lift the X-Y stage off the LED stand. Remove the top plate and its screws from the stage, and temporarily set aside. Lay the clean acetate template, illustrated in Figure 21, on the stage. The alignment box on the template must be carefully lined up with the rectangular cutout in the stage, which is outlined in blue in Figure 22. Tape the template in place, and fit the stage back onto the base of the LED stand. Gently lay the glass stage plate over the acetate to ensure it remains flat.

Using the eyepieces and moving the stage only along the x-axis, check the alignment of the template's horizontal ruler. If the ruler does not remain straight during the translation of the X-stage, make the necessary adjustments to the position of the template. Repeat for the y-axis using the vertical ruler. The stage will be reassembled after fine alignment.




Figure 20. Disassembling X-Y stage




Figure 21. PET acetate template drawing




Figure 22. Rough alignment



The aLIGO Particulate Evaluation Tool (PET)

5.3.6 Total Magnification

The total magnification of the microscope depends on the method of observation. For binocular observation (using the eyepieces), the total magnification is the product of the zoom magnification, ocular (eyepiece), and objective lenses:

$$M_{bino} = M_{obj} \times M_{oc} \times M_{zm}$$

M_{bino} : Total mag. for binocular observation

M_{obj} : Objective lens magnification

M_{oc} : Eyepiece magnification

M_{zm} : Zoom magnification

When using the digital camera, the eyepiece magnification is disregarded:

$$M_{cam} = M_{obj} \times M_{zm} \times M_{adpt}$$

M_{cam} : Total mag. for camera observation

M_{adpt} : Camera adapter magnification

For this setup, the total magnification will be the same whether the eyepieces or camera is used.

5.3.7 Microscope Settings Summary

The following table is a summary of the basic settings for the initial set up and alignment of the microscope. The same settings will also be used for swipe tool sample analysis.

Table 3. Basic AmScope settings


Setting	Value
Upper LED intensity	Max
Lower LED intensity	Max
Working distance	~ 4"
Objective lens	1X
Eyepieces	10X
Camera adapter	10X
Zoom magnification	1X
Total magnification – binocular	10X
Total magnification – camera	10X

The aLIGO Particulate Evaluation Tool (PET)

5.4 Software Setup

A software program called ToupView (64 bit version) is used to view samples through the digital camera and take particle count measurements. Before analyzing filter samples, ToupView must first be set up and calibrated for the swipe tool. If the software has not already been installed, do so now. Make sure the digital camera is connected to the computer's USB port. If any problems are encountered, refer to Section 9 for troubleshooting assistance.

5.4.1 ToupView Preferences

- 1) Open ToupView by double clicking on the  icon or locate ToupView in the Programs menu.
- 2) Choose **Options > Preferences** from the **Menu toolbar**. The **Preferences dialog** will appear, as shown in Figure 23.
- 3) Select the **Capture** page, and uncheck **Capture with Marker and Watermark**.
- 4) Select the **Misc** page, and check the option **Automatically activate the camera sidebar when startup**.
- 5) Select the **Rulers and Grids** page, and set **Style** to Auto Grid, **Line Style** to Solid, and **Color** to Gray.
- 6) Click **Ok**.




Figure 23. ToupView window and Preferences dialog

5.4.2 Camera Settings

The settings in this section are applicable only for the initial set up and calibration of the AmScope and ToupView software. Some of the settings will be adjusted for the analysis of actual swipe samples.

The aLIGO Particulate Evaluation Tool (PET)

- 1) Click **View > Sidebar > Camera**. There are 11 groups of settings in the **Camera Sidebar** which can be expanded by clicking on the group name.
- 2) Expand the **Camera List** group, and select **UCM50900KPB**. A live view from the microscope will appear, as shown in Figure 24.




Figure 24. ToupView live camera view and Camera Sidebar

- 3) Adjust the **Camera Sidebar** settings to match the values in Table 4.

Table 4. Camera Sidebar settings for calibration

Group	Setting	Value
Capture & Resolution	Live & Snap	3488*2616
Exposure & Gain	Auto Exposure	Checked
White Balance	Temperature Tint	6967 886
Color Adjustment	All settings	Default
Frame Rate	Frame Rate	40%
Color/Gray	Color/Gray	Color
Sampling	Sampling	Bin
Power Frequency	Power Frequency	AC (60 Hz)

The aLIGO Particulate Evaluation Tool (PET)

Note: If the video is still black, try decreasing the Frame Rate setting even further. Matching the [White Balance](#) values in [Table 1](#) can be difficult. The [White Balance](#) settings will be re-calibrated later, so only approximate values are required at this point.

- Expand the [Exposure & Gain](#) group, and check the [Auto Exposure](#) box. A green rectangle marking the reference region will appear over the video. Resize the rectangle to cover the entire video, then click [Default](#).

Note: The live view may momentarily go black, but should recalibrate itself within a few seconds. The setting values will vary depending on the lighting conditions in the room. As long as the setup remains consistent, the results should be repeatable.

- Depending on the orientation of the digital camera, it may be helpful to change the [Flip](#) settings. This is optional. The Caltech PET setup uses the [Vertical Flip](#) setting.

5.4.3 AmScope Fine Alignment

The live video feed can now be used to improve the alignment of the X-Y stage with the digital camera, focusing arm, and microscope body. The eyepieces and camera are slightly offset, so further adjustments to the AmScope setup will likely be required.

- Use the X-Y stage to bring the horizontal alignment ruler on the acetate template into the view in ToupView.
- Turn the x-axis knob on the stage, and watch the ruler in the live camera view as the stage travels along the entire range of the x-axis. The gridlines are a useful position reference.
- If the ruler is skewed in the live view, like [Figure 26](#), then make adjustments to the position of the focusing arm and microscope body as described in [Section 5.3.4](#).




Figure 25. Template ruler well aligned




Figure 26. Template ruler poorly aligned

- When the ruler remains perfectly straight along the entire x-axis, just as in [Figure 25](#), then these steps can be repeated to align the y-axis.
- Use the X-Y stage to find the filter alignment circle on the acetate template in ToupView. Both the top and bottom edges of the circle must be visible. If both edges cannot be seen, adjust the position of the template on the stage and return to the first step in this section.



The aLIGO Particulate Evaluation Tool (PET)

- 6) Remove the X-Y stage from the LED stand, taking care not to disturb the position of the acetate template.
- 7) Poke holes in the acetate to make room for the screws. Replace the stage's top plate, glass plate, plate stop, and screws.

5.4.4 Calibration

A magnification of 10X is used for the PET. If the microscope is needed for applications beyond swipe tool sample analysis, the steps in this section can be repeated for additional magnifications. Refer back to Section 5.3.6 for help calculating total magnification. If ToupView is not calibrated for the total magnification being used for the microscope, any measurements made may be inaccurate.

- 1) On the **ToupView toolbar**, adjust the **Unit** to Pixel, the **Magnification** to NA.
- 2) Go to **View > Best Fit**. This should change the **Zoom** setting to 32% if the ToupView screen is maximized.
- 3) Place the calibration slide, depicted in Figure 27, on the X-Y stage. Align the slide lengthwise on the x-axis.
- 4) Move the stage until the "1 DIV = 0.1 mm" scale is centered on screen, as demonstrated in Figure 28.
- 5) Bring the scale into focus using the focusing knob on the microscope.



Figure 27. Calibration slide




Figure 28. Calibration slide 0.1 mm scale in ToupView

The aLIGO Particulate Evaluation Tool (PET)

Note: Focus on the center of the scale. This will cause the ends to be slightly out of focus. Refer to Section 9.2 for help on focusing.

6) Change the **Zoom** to 100%, and center the 0.1 mm scale on the screen using the scrollbars in ToupView.

7) Click **Options** > **Calibrate** or click on the (Calibrate) icon in the ToupView toolbar. The **Calibrate dialog** and a red line with a number of pixels and 0.000 μm are overlaid on the video window. This can be seen in Figure 29.

8) Use the mouse to align the red line between two division marks on the calibration scale. Drag the ends of the red line as far as possible across the screen for the most accurate calibration, as in Figure 30.

9) In the **Calibrate dialog**, enter 10X into the **Magnification** field and the appropriate value in mm to the **Actual Length** field. In Figure 30, the **Actual Length** of the red line is 3 mm long. The value entered for the **Actual Length** will be displayed above the red line, and a value will be automatically calculated and displayed in the **Resolution** field.

Note: The resolution will be required later to convert the particle areas from pixels to μm .

10) Click **Ok**. ToupView is now calibrated for 10X magnification.

11) To view magnification and resolution settings, go to **Options** > **Magnification**.




Figure 29. Default Calibrate dialog




Figure 30. Calibration settings for 10X magnification



The aLIGO Particulate Evaluation Tool (PET)

6 SAMPLE ANALYSIS PROCEDURE

This section describes the process of obtaining a particle count from a filter sample. It is assumed the hardware and software have been properly calibrated. First time users should analyze the “dirty” control sample to practice the procedure before moving on to the other control samples or any samples collected on site. If any problems are encountered, refer to Section 9 for troubleshooting assistance.


6.1 Required Parts and Tools

In addition to the stereomicroscopy and image analysis system described in Section 5, the items listed in Table 5 are required for the analysis of swipe tool samples.

Table 5. Parts required for sample analysis

Part	P/N	Vendor	Purchasing Responsibility
Filter forceps	XX6200006P (3/pack)	Millipore , Amazon	<i>Sites to re-order as required</i>
PET custom microscope templates	LIGO-D1301007		
Acetate	3R3028 (100/pack)	Amazon	<i>Sites to re-order as required</i>
Printer	N/A	N/A	<i>Sites (general lab supplies)</i>
Guillotine paper cutter	N/A	N/A	<i>Sites (general lab supplies)</i>
Microscope slides	206A2 (72/pack)	Amazon	<i>Sites to re-order as required</i>
Vectra Alpha 10 wipes	TX1010	VWR	<i>Sites (general lab supplies)</i>
IPA	N/A	N/A	<i>Sites (general lab supplies)</i>

6.2 Sample Preparation and Setup

- 1) Wipe down the X-Y stage’s glass plate and a microscope slide with an IPA soaked Vectra Alpha 10 wipe, and allow both to dry.
- 2) Open ToupView by double clicking on the  icon or locate ToupView in the Program menu.
- 3) Select [UCMO50900KPB](#) in the [Camera Sidebar](#) to start the live video feed.
- 4) Make sure the current [Camera Sidebar](#) settings still match the calibration settings in Table 4.
- 5) Go to [View > Best Fit](#). This should change the [Zoom](#) setting to 32% if the ToupView screen is maximized.
- 6) On the [ToupView toolbar](#), set the [Unit](#) to Micrometer (μm) and [Magnification](#) to 10X.
- 7) Transfer the filter from the PetriSlide to X-Y stage using clean forceps, and center it over the filter alignment circle on the acetate template. Take care not to touch the sample area or disrupt the collected particulates during the transfer.

Note: Do not put a sample on the X-Y stage if it is still wet.

The aLIGO Particulate Evaluation Tool (PET)

- 8) Align the length of the filter's sample collection area with the stage's y-axis. As Figure 32 illustrates, samples typically have a bend line (or two) where the filter was folded over the swipe tool's sponge. Use the forceps to rotate the filter so the bend line is parallel to the stage's y-axis. The filter must still be centered over the circle on the acetate template.

Note: If the system was properly aligned in Section 5.4.3 and the filter is centered in the alignment circle, then both the top and bottom edges of a filter will always be visible in ToupView. It is good practice for first time users to confirm this before continuing.




Figure 31. Filter on X-Y stage centered over lower LED and covered with a microscope slide




Figure 32. Illustration of typical filter sample (not to scale)

- 9) Inspect the sample collection area for contamination. If there are numerous particles and fibers (~20 or more) distributed within the area, which is usually the case, skip to Step 11. However, if there are few or no particles, continue onto the next step.
- 10) Center a clean acetate stitching aid (from the PET custom microscope templates) over the sample collection area as per Figure 33.
- 11) Set the clean microscope slide directly over the filter (and stitching aid if used) to flatten and protect the sample area. The microscope slide should be roughly centered over the sample area, and parallel to the bend line.
- 12) Keeping the sample area in the live camera view, move the stage up or down the y-axis to find an area away from the edges of the filter that has multiple particulates. Bring the particles near the center of the live view into sharp focus.

Note: For cleaner samples, center one of the stitching aid's lines in the live view, and bring it into sharp focus.

- 13) Move the stage away from the sample collection area along the x-axis to locate a large un-swiped section of the filter, as demonstrated in Figure 34 on the following page. The selected area must be away from the edges of the filter, and must also be completely covered by the microscope slide.




Figure 33. Illustration of clean filter sample (not to scale)

The aLIGO Particulate Evaluation Tool (PET)

- 14) Go to [Camera Sidebar > White Balance](#). A red rectangle labelled White Balance is overlaid on the live view to mark the reference region.
- 15) Use the mouse to drag the red rectangle to cover as much blank filter area as possible. Click the [White Balance](#) button to recalibrate the color settings, and then collapse the group.
- 16) Move the stage to bring the sample collection area back into view. Select a location away from the edges of the filter.
- 17) Go to [Camera Sidebar > Exposure & Gain](#). A green rectangle labelled Exposure appears over the video to mark the reference region.
- 18) Expand the green rectangle to cover the entire live view. Click the [Default](#) button. The view may momentarily go black, but should recalibrate itself within a few seconds. Allow the [Exposure Target](#), [Exposure Time](#), and [Gain](#) values to stabilize.
- 19) Uncheck the [Auto Exposure](#) box, and collapse the group.




Figure 34. Recalibrating White Balance on un-swiped area




Figure 35. Recalibrating Exposure on sample collection area

6.3 Stitching

Multiple overlapping images of a sample are taken, and then combined into a single image. The stitching algorithm detects and matches features between images. Cleaner filters are a special case because there are insufficient common features in the overlapping fields. Without a stitching aid, images of samples with few or no particulates will not stitch successfully.

The stitching process can be challenging to get the hang of, and even harder to troubleshoot. The basic method should be attempted first, and if successful, the user may switch to the advanced method. Diagnosing causes of failure is easier in the basic method, but the process takes longer. It is important to note that if an error occurs during the advanced method, the entire stitching process simply fails and images taken of the sample cannot be recovered.

6.3.1 Basic Method

- 1) Click on the  ([Image select](#)) icon in the [ToupView toolbar](#). The cursor will turn into a crosshair.

The aLIGO Particulate Evaluation Tool (PET)

- 2) Draw a rectangle over the live image as shown in Figure 36. This rectangle should follow the gridlines on the screen; excluding 1 row of boxes on the top and on the bottom, 2 columns on the right, and 3 columns on the left.




Figure 36. Image select rectangle




Figure 37. Vignetting

Note: Drawing the rectangle is a way to work around vignetting, which has a huge impact on counting. Figure 37 is the same view as Figure 36, just with higher contrast to highlight vignetting. Drawing a rectangle around the least affected area will crop the worst sections when taking an image.

- 3) Re-focus on the particles at the center of the rectangle instead of the center of the live camera view.
- 4) Use the stage to move the live image up the y-axis until the top edge of the filter is just within the rectangle.
- 5) Allow the live view to stabilize, and click **Snap** in the **Camera Sidebar**. A tab for the cropped image, labelled 001, will appear beside the camera tab.




Figure 38. Preparing to Snap first image

Note: Blurry images are a very common source of stitching errors. When the X-Y stage is moved, the live view in ToupView may update in segments like an inch worm moving up the screen. Make sure the entire live view is stable before taking an image.

- 6) Rotate the Y translation knob counter-clockwise 10 increments to move the live view slightly down the sample. The graduation locator markings are clearly labelled on the knob itself. The dots on the X translation knob, depicted in Figure 39, are a useful reference.
- 7) Allow the picture to stabilize, and **Snap** another image.
- 8) Continue to move down the y-axis in increments of 10, and **Snap** images until the bottom edge of the filter comes into view.
- 9) As with the first image of the sample, move the stage down the y-axis until the bottom edge of the filter is within the rectangle, and **Snap**. The Y translation knob does not have to be moved the full 10 increments for the final image. There should be 14 images in total if the filter was centered on the swipe tool sponge.



Figure 39. X-Y stage knob

The aLIGO Particulate Evaluation Tool (PET)

- 10) Click on any of the image tabs. From the [Menu toolbar](#), select [Process](#) > [Stitch](#) to open the [Stitch dialog](#) shown in Figure 40.

Note: A different version of the Stitch dialog, described in the advanced method, will open if the live camera view is still open instead of one of the images.

- 11) In the [Stitch dialog](#), click [Add all](#) to select all of the photos. Choose [Vertical](#) for [Direction](#), [Foreground](#) for [Run](#), and click [Ok](#). A stitched image, rotated 90 degrees counter-clockwise, will appear in a new tab.

Note: This step can take a long time. Be patient.

- 12) Make sure the stitched image includes all of the images taken of the sample. The easiest way to do this is to look for the edges of the filter, as emphasized in Figure 41.

- 13) [Save](#) the image.

Note: Saving an image for the first time or selecting [Save As](#) will automatically close the image.




Figure 40. Stitch dialog in basic stitching method





Figure 41. Resulting stitched image with the 1st and 14th images highlighted

6.3.2 Advanced Method

- 1) Click on the  icon, and draw a rectangle over the live image to exclude 1 row of boxes on the top and bottom, 2 columns on the right, and 3 columns on the left
- 2) Use the stage to move the live image up the y-axis until the top edge of the filter is just barely within the rectangle.
- 3) With the live camera view selected, select [Process](#) > [Stitch](#) to open an alternate [Stitch dialog](#).

The aLIGO Particulate Evaluation Tool (PET)

- 4) Set the **Direction** to Vertical and **Run** to Foreground, as shown in Figure 42.
- 5) Click **Snap** to capture the current image. The image will appear in the **Stitch dialog**.
- 6) Rotate the y-axis knob counter-clockwise 10 increments, allow the picture to stabilize, and **Snap** another image.
- 7) Continue to move down the y-axis in increments of 10, and **Snap** images until the bottom edge of the filter comes into view.
- 8) Move the stage down the y-axis until the bottom edge of the filter is just within the rectangle. The stage does not have to be moved the full 10 increments for the final image.
- 9) Click **Ok**.

Note: This step can take a long time. Be patient.

- 10) Make sure the stitched image includes all of the images taken of the sample. The easiest way to do this is to look for the edges of the filter.
- 11) **Save** the image.

6.4 Image Processing

- 1) Go to **File > Open Image** and select the stitched image that was just saved.
 - 2) Go to **File > Save As** to make a copy of the original image. For every sample, an unprocessed and processed stitched image must be saved.
- Note: Saving the original images will allow the PET analysts to re-analyze samples without starting from scratch should the image processing procedure evolve over time.*
- 3) In the **Menu toolbar**, select **Process > Filter** to open the **Filters dialog**.

- 4) Under the **Image Enhance** tab; set the **Filter** to High Gauss, **Option** to 7 x 7, **Passes** to 1, and **Strength** to 1.
- 5) Click **Ok** to apply the filter.

Note: It is helpful to have the original image available for reference for this and the remaining steps in this section. Open the image in a different program, such as Windows Photo Viewer.




Figure 42. Stitch dialog in advanced stitching method





Figure 43. Filters dialog with High Gauss settings

The aLIGO Particulate Evaluation Tool (PET)



Figure 44. Drawing crop box around stitching image

- 6) Click on the  icon, and draw a rectangle over as much of the stitched image as possible, excluding any white areas from the LEDs and black edges around images. Zoom in if necessary.
- 7) In the **Menu toolbar**, select **Image > Crop**.
- 8) Go to **Image > Mode > Gray Scale**.

9) Go to **Process > Segmentation** to open the **Segmentation dialog**. This thresholding technique is used to separate darker particles from the filter background.




Figure 45. Segmentation (red) with an excessively large threshold value

10) Drag the left threshold bar in the **Segmentation dialog** all the way left to 0. Pixels that are true black, i.e. equal to 0, are highlighted in a different color. The highlight color is red in Figure 45.




Figure 46. Segmentation (white) correctly highlighting dark particles

11) Raising the value of the right threshold bar will highlight progressively lighter colored pixels. Adjust the right threshold bar so that only the dark particles are highlighted. The threshold value is too high if highlighted pixels appear where there were no particles in the original stitched image, like shadowing around the edges of the image in Figure 45. Figure 46 is an example of how to correctly apply the segmentation process. The optimal threshold range will vary depending on the microscope setup.

- 12) Switch the highlight color from red to white using the drop down box in the upper left corner of the **Segmentation dialog**, and click **Ok**.
- 13) Select **Process > Filter** to open the Filters dialog again.

The aLIGO Particulate Evaluation Tool (PET)

- 14) Under the **Morphological** tab; set the **Filter** to Open, **Option** to 5 x 5 Circle, and **Passes** to 1. For images with bright objects on a darker field, this filter will smooth object contours and help fill in small holes.

Note: Segmentation, while extremely useful, can inadvertently modify the shape of objects. For example, it may erode the thinnest point of a fiber, making it appear as though there are two smaller fibers in the final image. The open filter helps to minimize this potential source of error.

- 15) Click **Ok**.
- 16) Go to **Process > Binary** to open the **Binary dialog**.
- 17) Drag the binary threshold bar toward the right or use the arrow keys. If the value of a pixel is greater than the selected threshold, the pixel will be turned white. Otherwise, it will be turned black. Adjust the threshold value so that the background is completely black and the particles are white. If the threshold is moved too far to the right, the particles will begin to shrink or disappear. The optimal threshold value will vary by setup. The maximum threshold value, 255, was used for the sample in Figure 48.
- 18) Click **Ok**.
- 19) **Save** the processed image.




Figure 47. Filters dialog with Open filter settings






Figure 48. Conversion to a binary image

The aLIGO Particulate Evaluation Tool (PET)

6.5 Image Size

- 1) Go to **Image > Image Scale** to open the **Scale Image dialog** with the height and width in pixels. Make a note of the dimensions, as they will be required later.

Note: First time users should complete the remaining steps in this section to confirm ToupView is properly calibrated. Experienced users should skip to Section 6.6.

- 2) Select the  (**Line**) icon from the **ToupView toolbar**, or go to **Measurements > Line > Arbitrary Line**.
- 3) Draw a straight vertical line from the top of the stitched image to the bottom, as in Figure 49. The length in μm of the line will appear in the **Measurement Sidebar**.
- 4) Draw a straight horizontal line to measure the width of the image.
- 5) Confirm the image dimensions in the **Scale Image dialog** match the user-measured dimensions. This can be done by dividing the **Scale Image dialog** dimensions (pixels) by the resolution (pixels/meter). Alternatively, the **Unit** can be changed from Micrometer (μm) to Pixel in the **ToupView toolbar** to automatically convert the user-measured length.
- 6) To remove the measurement lines, click on the  (**Object Select**) icon or go to **Measurements > Object Select**.
- 7) Select the red lines by clicking on them, and then go to **Edit > Cut**.




Figure 49. Verifying image dimensions

The aLIGO Particulate Evaluation Tool (PET)

6.6 Particle Count

- 1) Select **Plugin > Count** from the **ToupView toolbar**. A new **Count Window** will be open.
*Note: The **Count Window** will open with whatever image is currently selected or in view in the main **ToupView Window**. Only one image can be opened in the **Count Window** at a time.*
- 2) Select **Settings > Parameter** to open the **Setting dialog**.
- 3) Enter the settings defined in Figure 50 and Table 6. The user may select any color for the **Border**, **Number**, and **Center** fields. The gray min and max values determine which pixels will be counted. Since black is equal to 0, the software will count everything else.




Figure 50. ToupView Count Window and Setting dialog

Table 6. Count settings

Setting	Value
GrayMin	1
GrayMax	256
AreaMin	1
AreaMax	1,000,000,000
PeriMin	1
PeriMax	1,000,000,000

- 4) Click **Ok**.
- 5) Go to **M&C > Count**. Each particle will automatically be outlined and labelled with a unique number.

The aLIGO Particulate Evaluation Tool (PET)



Figure 52. Section of binary image before count




Figure 51. Section of binary image after count

- 6) Use the icon to zoom in, and look closely at the particles. Verify the highlighted outlines match the true perimeters of the particles.
- 7) Go to **M&C > Check the Result...** to open the **Results dialog**. The counted particles are listed by number with their corresponding area and perimeter values.
- 8) Click **Export** to open the count data in Excel. The area and perimeter values are listed in Column A and Column B respectively. The particle number (SN) is not exported.
- 9) **Save** the Excel spreadsheet.




Figure 53. Count results dialog



The aLIGO Particulate Evaluation Tool (PET)

6.7 Summary of Steps and Settings

Whenever a new sample is analyzed, the Camera Sidebar settings should be reset to the values that were used to calibrate ToupView. The Exposure & Gain and the White Balance group setting are then calibrated for each individual sample. Table 7 and Table 8 summarize the software settings used for image analysis and counting.

Figure 54 on the subsequent page illustrates the main steps in the image analysis procedure.

Table 7. Initial and final Camera Sidebar settings for sample analysis

Group	Setting	Initial Value	Final Value
Capture & Resolution	Live & Snap	3488*2616	3488*2616
Exposure & Gain	Auto Exposure Exposure Target & Time Gain	Checked (Auto) (Auto)	<i>Not Checked</i> <i>Recalibrated</i> <i>Recalibrated</i>
White Balance	Temperature Tint	6967 886	<i>Recalibrated</i> <i>Recalibrated</i>
Color Adjustment	All settings	Default	Default
Frame Rate	Frame Rate	40%	40%
Color/Gray	Color/Gray	Color	Color
Sampling	Sampling	Bin	Bin
Power Frequency	Power Frequency	AC (60 Hz)	AC (60 Hz)

Table 8. Image processing and counting settings

Process	Setting	Value
High Gauss filter	Option (kernel)	7 x 7
	Passes	1
	Strength	1
Segmentation	Upper threshold	0
	Lower threshold	Determined by user
Open filter	Option (kernel)	5 x 5
	Passes	1
Binary	Threshold	Determined by user
Count	GrayMin	1
	GrayMax	256
	AreaMin	1
	AreaMax	1,000,000,000
	PeriMin	1
	PeriMax	1,000,000,000



The aLIGO Particulate Evaluation Tool (PET)




Figure 54. Summary of steps in image analysis procedure



The aLIGO Particulate Evaluation Tool (PET)

7 PARTICLE CLEANLINESS LEVEL CALCULATOR

The final, and easiest, step in the PET System is to calculate the PCL using an Excel template, [LIGO-E1300641](#). This section gives instructions for using the spreadsheet, and explains how the PCL is calculated.

100 samples can be entered into the PCL calculator template. For each observatory, a separate spreadsheet will be needed for every chamber and main work area, such as the LVEA, LBR, etc. Similar to the FBI samples and 4” wafers, all swipe tool sample data will be maintained on the DCC in a top-level library, [LIGO-T1300921](#).

Figure 55 shows the tabs for the PCL calculator’s eight sheets. If sample data is entered correctly into the Raw Data sheet, the template will automatically calculate and graph PCLs. *Do not make any changes to the template’s other sheets.*




Figure 55. Excel PCL calculator tabs

7.1 Raw Data Entry & Results

- 1) Check the top-level library, [LIGO-T1300921](#), to see if there is already a spreadsheet dedicated to the area where the sample was collected. If so, download that spreadsheet from the DCC, and skip to Step 3. If not, download the latest version of the template, [LIGO-E1300641](#).
- 2) To set up a new spreadsheet, enter the following into the Raw Data sheet:
 - Site
 - Location
 - Resolution (pixel/m)
 - Template version

Note: To look up the resolution in ToupView, go to *Options* > *Magnification* in the *ToupView* toolbar.

Raw Data							
SITE:	Caltech	SAMPLES					
LOCATION:	Contamination Lab	Sample Name	Example #1	Example #2	Example #3	Example #4	Example #5
RESOLUTION (pixel/m):	296267	Date	11/25/13	12/9/13	12/13/13	1/7/14	1/9/14
TEMPLATE VERSION:	v3	Swipe Length (in)	18	24	20	15	12
		Image Width (pixels)	13762	13750	13743	13776	13748
		Image Height (pixels)	2041	2000	2018	2036	2028
		Notes	Optics table, left side	Optics table, right side	Table near air knife experiment	Optics table, right side	Table near air knife experiment
		ToupView Raw Data - Area (pixels)	33	41	45	20	826
			155	50	200	65	1259
			604	32	631	1218	2486
			32	1	23	15	50
			1	69	61	1463	76
			478	168	54	440	93
			46	218	94	3249	60
			1	88	1	3	156
			1	71	1	1	1279
			1	48	5	1	101
			1	1	3	1	251
			1	19	34	1	2
			92	12	75	1	27

Figure 56. Raw Data sheet with setup data highlighted in red and ToupView counts highlighted in blue



The aLIGO Particulate Evaluation Tool (PET)

3) For each sample, manually enter the following into the Raw Data sheet:

- Sample name
- Date sample was collected
- Swipe length (inches)
- Stitched image height (pixels) after cropping
- Stitched image width (pixels) after cropping
- Notes (such as more detailed information on the location or any additional useful information)

*Note: To look up the image dimensions in ToupView, go to **Image > Image Scale** in the **ToupView toolbar**.*

4) For each sample, copy the particle areas (pixels) from Column A of the exported ToupView count data, and paste it into the Raw Data sheet.

5) Go to the Results sheet. The PCL, total number of particles, minimum particle diameter (μm), maximum diameter (μm), swipe sample area (m^2), and image area (m^2) are listed for each sample. There is also a graph tracking the changes in PCL over time.




Figure 57. Results sheet

6) Go to the appropriate graph sheet. The samples' PCLs are plotted on log-log² axes along with the standard PCLs (grey) and goal PCL (red), as shown in Figure 58 on the subsequent page.



The aLIGO Particulate Evaluation Tool (PET)




Figure 58. IEST-STD-CC1246D standard PCLs and swipe tool sample PCLs plotted on log-log² axes

7) Save and upload the spreadsheet to the DCC. Make sure it is linked to the top-level swipe tool sample library.

At this point, the user will have completed all the required steps in the PET. The subsequent sections explain how the calculations are performed.

7.2 Definition of Cleanliness Levels

IEST-STD-CC1246D defines a uniform set of particle size distribution criteria for surfaces with the following equation:

$$\log N = 0.926(\log^2 L - \log^2 x)$$

N : number of particles per 0.1 m² of surface

0.926: absolute value of the slope of the line

L : cleanliness level (aka PCL)

x : particle diameter (µm)



The aLIGO Particulate Evaluation Tool (PET)

Levels are named for the largest particle allowed per specified area. For example, a single 30 μm particle on a 0.1m² surface corresponds to PCL 30.

7.3 Calculations

There is an underlined title above each table on the Calculations sheet explaining the purpose of the computation performed, as shown in Figure 59.

7.3.1 Areas

The swipe area and image area are calculated from the user entered data as follows:

$$A_{swipe} = [l(\text{in}) \times 0.0254(\text{m/in})] \times \left[\frac{w(\text{pixels})}{R(\text{pixels/m})} \right]$$

$$A_{image} = \frac{w(\text{pixels}) \times h(\text{pixels})}{[R(\text{pixels/m})]^2}$$

A_{swipe} : swipe area

A_{image} : image area

l : swipe length

w : stitched image width after cropping

R : microscope resolution for 10X magnification

h : stitched image height after cropping

	A	C	D	E	F	G	H
1	Calculations						
2							
3							
4							
5	<u>Calculates areas:</u>						
6			Example #1	Example #2	Example #3	Example #4	Example #5
7		Swipe area (m ²)	2.12E-02	2.83E-02	2.36E-02	1.77E-02	1.41E-02
8		Image area (m ²)	3.20E-04	3.13E-04	3.16E-04	3.20E-04	3.18E-04
9							
10	<u>Bins data:</u>						
11		Particles per 0.1 m ² , N					
12		Particle diameter (μm), x	Example #1	Example #2	Example #3	Example #4	Example #5
13		5	9,808	1,909	59	6,000	12,747
14		10	7,167	1,145	51	4,916	10,669
15		15	5,118	884	51	4,019	8,802
16		20	3,682	707	47	3,280	6,985
17		25	2,609	573	42	2,726	5,508
18		30	1,808	449	30	2,230	4,355
19		35	1,356	364	21	1,868	3,429
20		40	998	269	17	1,513	2,786
21		45	782	216	17	1,236	2,270
22		50	593	180	13	1,010	1,930

Figure 59. Area and binning sections on Calculations sheet



The aLIGO Particulate Evaluation Tool (PET)

7.3.2 Binning

The second table on the Calculations, lists particle diameters from 5 to 500 μm in increments of 5 μm. For a given diameter (μm), the corresponding area (pixels) is calculated assuming circular particles.

$$A_p = \pi[0.5x(\text{pixels}) \times R(\text{pixels/m})]^2$$

A_p : particle area

For a given diameter (μm), the table counts and bins every particle from the ToupView raw data greater than or equal to the corresponding area (pixels).

7.3.3 PCL

Using the IEST-STD-CC1246D equation from Section 7.2, the PCL is calculated for each particle diameter (μm). The largest PCL is the overall PCL for the sample. The cells for the largest PCL values are colored red to highlight the particle diameters responsible for the overall PCL.

Calculates PCL:		PCL, L				
Particle diameter (μm), x	Example #1	Example #2	Example #3	Example #4	Example #5	
5	155	102	35	137	165	
10	187	119	49	171	206	
15	210	137	63	198	237	
20	226	153	74	220	261	
25	237	166	84	240	281	
30	244	176	88	256	297	
35	252	186	90	271	309	
40	257	190	94	282	322	
45	263	196	104	291	332	
50	266	203	105	299	344	
55	271	213	100	310	355	
60	280	218	108	321	359	
65	286	213	116	327	367	
70	289	220	124	336	376	
75	294	223	132	344	378	
80	297	226	139	347	385	
85	301	235	147	360	390	
90	311	242	155	355	392	
95	312	240	137	364	397	
100	314	240		374	397	
105	318	238		385	405	
110	327	241		387	414	
115	335	241		390	424	
120	343	242		403	439	
125	351	250		412	444	
130	344	232		418	447	
135	355	240		422	458	
140	360	248		428	456	
145	342	256		437	465	
150	353	209		442	475	
155	355	216		450	475	
160	365	222		457	488	
165	364	229		469	495	
170	360	235		475	492	
175	370	242		487	492	
180	379	248		488	491	
185	376	255		493	495	
190	385	261		483	498	
195	394	268		487	491	
200	403	274		498	492	
205	411	281		508	476	
210	420	287		519	466	
215	429	294		509	496	
220	438			510	491	
225	447			507	461	
230	456			506	470	
235	443			515	452	
240	451			524	461	
245	460			534	469	
250	402			543	440	
255	409			535	388	
260	417			524	395	
265	424			513	402	
270	431			521	410	
275	439			530	417	
280	395			465	424	
285	402			472		
290	408			480		
295	415			487		
300	422			496		
305	428			503		
310	435			510		
315	442			518		
320	448			525		
325	455			472		
330	461			479		
335	468			486		
340						
345						
350						
355						
360						
365						
370						
375						
380						
385						
390						
395						
400						
405						
410						
415						
420						
425						
430						
435						
440						
445						
450						
455						
460						
465						
470						
475						
480						
485						
490						
495						
500						
Overall PCL	468	294	155	543	498	

Figure 60. PCL section on Calculations sheet



The aLIGO Particulate Evaluation Tool (PET)

8 SUMMARY OF SETTINGS

Settings described in the previous sections are repeated here for easy referencing.

8.1 Microscope

The same basic settings are used for the initial setup and alignment of the microscope, ToupView calibration, and swipe sample analysis.

Table 9. Basic AmScope settings

Setting	Value
Light intensity	Max
Working distance	4"
Objective lens	1X
Eyepieces	10X
Camera adapter	10X
Zoom magnification	10X
Total magnification – binocular	10X
Total magnification – camera	10X

8.2 ToupView

The initial Camera Sidebar settings listed in Table 10 are used to calibrate the microscope. Whenever a new sample is analyzed, the Camera Sidebar settings are reset to the initial values. The Exposure & Gain and the White Balance group setting are then calibrated for each individual sample.

Table 10. Initial and final Camera Sidebar settings for sample analysis

Group	Setting	Initial Value	Final Value
Capture & Resolution	Live & Snap	3488*2616	3488*2616
Exposure & Gain	Auto Exposure Exposure Target & Time Gain	Checked (Auto) (Auto)	<i>Not Checked</i> <i>(Recalibrated)</i> <i>(Recalibrated)</i>
White Balance	Temperature Tint	6967 886	<i>Recalibrated</i> <i>Recalibrated</i>
Color Adjustment	All settings	Default	Default
Frame Rate	Frame Rate	40%	40%
Color/Gray	Color/Gray	Color	Color
Sampling	Sampling	Bin	Bin
Power Frequency	Power Frequency	AC (60 Hz)	AC (60 Hz)



The aLIGO Particulate Evaluation Tool (PET)

Much of the image processing settings are standardized between the sites and Caltech. However, the segmentation and binary threshold values will be determined by the PET analysts at each location.

Table 11. Image processing and counting settings

Process	Setting	Value
High Gauss filter	Option (kernel)	7 x 7
	Passes	1
	Strength	1
Segmentation	Upper threshold	0
	Lower threshold	Determined by user
Open filter	Option (kernel)	5 x 5
	Passes	1
Binary	Threshold	Determined by user
Count	GrayMin	1
	GrayMax	256
	AreaMin	1
	AreaMax	1,000,000,000
	PeriMin	1
	PeriMax	1,000,000,000

The aLIGO Particulate Evaluation Tool (PET)

9 TROUBLESHOOTING

This section provides tips for solving common ToupView issues. A fix for many ToupView problems is restarting the program and, if that fails, restarting the computer.

9.1 Black Camera View

If the camera view is black, make sure the digital camera is connected to the computer's USB port, and the Light Path Lever on the microscope is set to "IN." Try reducing the [Frame Rate](#) in the [Camera Sidebar](#).

9.2 Focusing

Subtle changes in the focus can drastically affect the image, especially in the x-axis (this may be unique to the Caltech setup due to the orientation of the camera, uneven surface, etc.). Figure 61 is a series of images taken at 10X magnification of the same sample, in the same spot; only the focus was adjusted. The red boxes highlight the area that is currently in focus. Users should try to focus on particles in the center of the live view.




Figure 61. Shifting areas of focus on the same sample

The aLIGO Particulate Evaluation Tool (PET)

9.3 Incomplete Stitching

Figure 62 is an example of incomplete stitching, which is usually because of an issue with one or more images. Look for identifying markers on the partially stitched image, and identify which image(s) failed to stitch. Common causes are:

- Blurriness in a single image, as shown in Figure 63.
- Poor focus
- Auto Exposure was checked, which can change the lighting and contrast between the individual photos.
- Misaligned X-Y staged evidenced by black areas in the photos, as in Figure 64.




Figure 62. Color bands show changes in settings between images for this incomplete stitch



Figure 63. Blurry image




Figure 64. Image from misaligned stage



The aLIGO Particulate Evaluation Tool (PET)

9.4 Failed Stitch

Try stitching just two images to narrow down the cause of failure. Alternatively, use the Stitch box on the acetate X-Y stage template.

9.5 Count Plugin Will Not Open

Only one image can be opened at a time in the [Count Window](#). If the plugin is not already open, make sure there is an image open in the ToupView

9.6 Particles Lost/Gained

Adjust the segmentation and binary threshold values.

9.7 Toolbars Disappeared

If the Menu toolbar or ToupView toolbar are deleted or have disappeared, the entire program must be reset. Go to [Start](#) > [Computer](#) > [Local Disc \(C:\)](#) > [Program Files](#) > [AmScope](#) > [ToupView](#) > [x64](#) > [clean](#). The software will then have to be recalibrated.