

LIGO-E1101109-v1

PR303 Mass measurement at Caltech

November 17, 2011

GariLynn Billingsley

PR303 on fixture	19.812 Kg
Fixture alone	<u>-7.667 Kg</u>
PR303 Mass	12.145 Kg

For measurement accuracy please refer to the excel file found at LIGO-T1100174-v4

Look for the tab:

A&D GP 60K Page 83

Our vendor uses the scale GP 60K, absolute accuracy after calibration by internal mass (as found in the above file) is 5g.